


Anti-Prrx2 Rabbit Polyclonal Antibody For Wb, Elisa - Q99811

Item Number: CM0021457



Introduction

PRRX2 Rabbit pAb is a polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, rat PRRX2. Ideal for developmental biology and wound healing research.

[Learn More](#)

Parameter	Value
Product Name	PRRX2 Rabbit pAb
Remarks/Alias	PMX2; PRX2; PRRX2
Species	Human
GeneID (Human)	51450
GeneID	51450
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-253 of human PRRX2 (NP_057391.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99811
Protein Weight	27kDa
Shipping	Ice bag