

Anti-Prss12 Rabbit Polyclonal Antibody For Wb, Elisa - P56730

Item Number: CM0014424

Introduction

Rabbit polyclonal antibody targeting human PRSS12 (Neurotrypsin), validated for Western blot and ELISA. Shows cross-reactivity with mouse and rat. Suitable for neuronal plasticity and memory research.

[Learn More](#)

Parameter	Value
Product Name	PRSS12 Rabbit pAb
Remarks/Alias	MRT1; BSSP3; BSSP-3; PRSS12
Species	Human
GeneID (Human)	8492
GeneID	8492
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-300 of human PRSS12 (NP_003610.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P56730
Protein Weight	97kDa
Shipping	Ice bag