


Anti-Prx Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bxm0

Item Number: CM0020690



Introduction

Rabbit polyclonal antibody against human PRX, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience, myelin biology, and Charcot-Marie-Tooth disease research.

[Learn More](#)

Parameter	Value
Product Name	PRX Rabbit pAb
Remarks/Alias	CMT4F; PRX
Species	Human
Gene ID (Human)	57716
Gene ID	57716
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1250-1350 of human PRX (NP_870998.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXM0
Protein Weight	155kDa
Shipping	Ice bag