


Anti-Psa/Klk3 Rabbit Monoclonal Antibody For Wb, Elisa - P07288

Item Number: CM0009117



Introduction

Anti-PSA/KLK3 rabbit monoclonal antibody (clone A4355) validated for WB and ELISA. Targets human prostate-specific antigen (UniProt P07288, 29 kDa). Immunogen: synthetic peptide from aa 162-261. Suitable for cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	PSA/KLK3 Rabbit mAb
Remarks/Alias	APS; PSA; hK3; KLK2A1; PSA/KLK3
Species	Human
GeneID (human)	354
GeneID	354
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 162-261 of human PSA/KLK3 (P07288).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P07288
Protein Weight	29kDa
Shipping	Ice bag