


Anti-Psap (Prosaposin) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07602

Item Number: CM0022592



Introduction

Rabbit polyclonal antibody against human PSAP (Prosaposin), validated for WB, IHC-P, IF/ICC, ELISA. Detects a ~58 kDa band. Cross-reacts with mouse and rat. Ideal for neurodegeneration and lysosomal storage disease research.

[Learn More](#)

Parameter	Value
Product Name	PSAP Rabbit pAb
Remarks/Alias	GLBA; SAP1; SAP2; PSAPD; PARK24; PSAP
Species	Human
GeneID (Human)	5660
GeneID	5660
Immunogen	Recombinant protein corresponding to amino acids 275-524 of human PSAP (NP_002769.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07602
Protein Weight	58kDa
Shipping	Ice bag