


Anti-Psap Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, Elisa - P07602

Item Number: CM0006353



Introduction

High-quality rabbit monoclonal antibody against human PSAP (prosaposin), validated for WB, IHC-P, ELISA. Targets lysosomal sphingolipid activator protein involved in glycosphingolipid metabolism and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] PSAP Rabbit mAb
Remarks/Alias	GLBA; SAP1; SAP2; PSAPD; PARK24; [KD Validated] PSAP
Species	Human
Gene ID (Human)	5660
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 275-524 of human PSAP (NP_002769.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P07602
Protein Weight	58kDa
Shipping	Ice bag