


Anti-Psca Rabbit Polyclonal Antibody For Wb And Elisa - O43653

Item Number: CM0029011



Introduction

Rabbit polyclonal antibody targeting human Prostate Stem Cell Antigen (PSCA), suitable for Western blot and ELISA. Useful in cancer and Alzheimer disease research. Reliable IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	PSCA Rabbit pAb
Remarks/Alias	PRO232; IncPSCA; PSCA
Species	Human
GeneID (Human)	8000
GeneID	8000
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43653
Protein Weight	12kDa
Shipping	Ice bag