

Anti-Psd2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bqi7

Item Number: CM0015721

Introduction

Anti-PSD2 rabbit polyclonal antibody validated for WB and ELISA in human, mouse, and rat. Ideal for studying PH and SEC7 domain-containing protein 2's role in membrane trafficking and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	PSD2 Rabbit pAb
Remarks/Alias	EFA6C; PSD2
Species	Human
Gene ID (Human)	84249
Gene ID	84249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human PSD2 (NP_115665.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQI7
Protein Weight	85kDa
Shipping	Ice bag