

Anti-Psd93/Chapsyn-110/Dlg2 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q15700

Item Number: CM0002841



Introduction

Rabbit monoclonal antibody targeting PSD93/DLG2, suitable for WB, IF-P, and ELISA. Recognizes human, mouse, and rat PSD93/DLG2. Useful for research on NMDA receptor signaling, chronic pain, and synaptic stability.

[Learn More](#)

Parameter	Value
Product Name	PSD93/chapsyn-110/DLG2 Rabbit mAb
Remarks/Alias	PSD93; PSD-93; PPP1R58; chapsyn-110; PSD93/chapsyn-110/DLG2
Species	Human
Gene ID	1740
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of human PSD93/chapsyn-110/DLG2 (Q15700).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15700
Protein Weight	98kDa
Shipping	Ice bag