

Anti-Psd95 Polyclonal Antibody For Wb, If/Icc, Elisa - P78352

Item Number: CM0027285



Introduction

High-quality rabbit polyclonal antibody targeting human PSD95 (DLG4, SAP90). Validated for WB, IF/ICC, ELISA with human, mouse, rat cross-reactivity. Ideal for synaptic plasticity and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PSD95 Rabbit pAb
Remarks/Alias	MRD62; PSD95; SAP90; SAP-90
Species	Human
GeneID (Human)	1742
GeneID	1742
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 420-510 of human PSD95 (NP_001308004.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78352
Protein Weight	80 kDa
Shipping	Ice bag