

Anti-Psd95 Rabbit Polyclonal Antibody For Wb, Elisa - P78352

Item Number: CM0028665



Introduction

Rabbit polyclonal antibody targeting human PSD95 (Disks large homolog 4, SAP90). Suitable for Western Blot and ELISA applications; cross-reacts with human, mouse, and rat samples. Essential for studying synaptic plasticity and postsynaptic density organization.

[Learn More](#)

Parameter	Value
Product Name	PSD95 Rabbit pAb
Remarks/Alias	MRD62; PSD95; SAP90; SAP-90
Species	Human
GeneID (Human)	1742
GeneID	1742
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PSD95 (NP_001308004.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78352
Protein Weight	80kDa
Shipping	Ice bag