

Anti-Psd95 Rabbit Polymab® For Wb, Ihc-P, Elisa - P78352

Item Number: CM0012016



Introduction

High-quality rabbit monoclonal PSD95 antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for synaptic plasticity and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PSD95 Rabbit PolymAb®
Remarks/Alias	MRD62; PSD95; SAP90; SAP-90
Species	Human
GeneID (Human)	1742
GeneID	1742
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PSD95 (NP_001308004.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P78352
Protein Weight	80kDa
Shipping	Ice bag