

Anti-Psen2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P49810

Item Number: CM0026685



Introduction

High-quality rabbit polyclonal antibody against human PSEN2, suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat PSEN2, ideal for Alzheimer's and gamma-secretase research.

[Learn More](#)

Parameter	Value
Product Name	PSEN2 Rabbit pAb
Remarks/Alias	AD4; PS2; AD3L; STM2; CMD1V; PSEN2
Species	Human
GeneID (human)	5664
GeneID	5664
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 267-384 of human PSEN2 (NP_000438.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49810
Protein Weight	50kDa
Shipping	Ice bag