

Anti-Psg1 Polyclonal Antibody For Wb/If-P/Ihc-P/Elisa - P11464

Item Number: CM0027965

Introduction

Rabbit polyclonal antibody against human PSG1 (pregnancy-specific beta-1-glycoprotein 1). Suitable for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects a ~47 kDa secreted glycoprotein.

[Learn More](#)

Parameter	Value
Product Name	PSG1 Rabbit pAb
Remarks/Alias	SP1; B1G1; PBG1; CD66f; PSBG1; PSG95; PSGGA; DHFRP2; PSBG-1; PSGIIA; FL-NCA-1/2; PS-beta-C/D; PS-beta-G-1; PSG1
Species	Human
GeneID (Human)	5669
GeneID	5669
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-419 of human PSG1 (NP_001171754.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P11464
Protein Weight	47kDa
Shipping	Ice bag