


Anti-Psip1 Rabbit Polyclonal Antibody For Wb, Elisa - 075475

Item Number: CM0019047



Introduction

Rabbit polyclonal antibody against human PSIP1 (LEDGF/p75), validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for research in transcriptional regulation, lens development, and viral integration.

[Learn More](#)

Parameter	Value
Product Name	PSIP1 Rabbit pAb
Remarks/Alias	p52; p75; PAIP; DFS70; LEDGF; PSIP2; PSIP1
Species	Human
Gene ID (Human)	11168
Gene ID	11168
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-530 of human PSIP1 (NP_150091.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75475
Protein Weight	60kDa
Shipping	Ice bag