


Anti-Psma1 Polyclonal Antibody For Wb - P25786

Item Number: CM0021700



Introduction

High-quality rabbit polyclonal antibody against PSMA1 (Proteasome subunit alpha type-1), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for proteasome function and protein degradation studies.

[Learn More](#)

Parameter	Value
Product Name	PSMA1 Rabbit pAb
Remarks/Alias	NU; HC2; PROS30; HEL-S-275; PSMA1
Species	Human
GeneID (Human)	5682
GeneID	5682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-263 of human PSMA1 (NP_002777.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25786
Protein Weight	30kDa
Shipping	Ice bag