


Anti-Psma2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P25787

Item Number: CM0004584



Introduction

Rabbit monoclonal antibody against human proteasome subunit alpha type-2 (PSMA2), suitable for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous PSMA2 in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	PSMA2 Rabbit mAb
Remarks/Alias	MU; HC3; PSC2; PMSA2; PSMA2
Species	Human
Gene ID (Human)	5683
Gene ID	5683
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PSMA2 (P25787)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25787
Protein Weight	26kDa
Shipping	Ice bag