


Anti-Psma2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P25787

Item Number: CM0021357



Introduction

High-quality PSMA2 Rabbit Polyclonal Antibody (P25787) validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat PSMA2. Essential tool for proteasome function and protein degradation studies.

[Learn More](#)

Parameter	Value
Product Name	PSMA2 Rabbit pAb
Remarks/Alias	MU; HC3; PSC2; PMSA2; PSMA2
Species	Human
Gene ID (Human)	5683
Gene ID	5683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human PSMA2 (NP_002778.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25787
Protein Weight	26kDa
Shipping	Ice bag