

Anti-Psma3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P25788

Item Number: CM0012362



Introduction

Anti-PSMA3 rabbit monoclonal antibody suitable for WB, IHC-P, ELISA with human, mouse, rat. Detects proteasome subunit alpha-3 (~28kDa), also known as C8, PSC3, HC8.

[Learn More](#)

Parameter	Value
Product Name	PSMA3 Rabbit mAb
Remarks/Alias	HC8; PSC3; PSMA3
Species	Human
GeneID (Human)	5684
GeneID	5684
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25788
Protein Weight	28kDa
Shipping	Ice bag