

Anti-Psma3 Rabbit Polyclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P25788

Item Number: CM0018132



Introduction

High-quality rabbit polyclonal antibody against human PSMA3, a 20S proteasome subunit vital for ubiquitin-dependent and independent protein degradation. Validated for WB, IP, IF-P, IHC-P, ELISA in human, mouse, rat. Ideal for proteasome research.

[Learn More](#)

Parameter	Value
Product Name	PSMA3 Rabbit pAb
Remarks/Alias	HC8; PSC3; PSMA3
Species	Human
Gene ID (Human)	5684
Gene ID	5684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human PSMA3 (NP_002779.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25788
Protein Weight	28kDa
Shipping	Ice bag