


Anti-Psma5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28066

Item Number: CM0028874



Introduction

Rabbit polyclonal antibody targeting human PSMA5 (proteasome subunit alpha type-5). Validated for WB, IHC-P, IF/ICC, ELISA with reactivity in human, mouse, rat. Ideal for proteasome research.

[Learn More](#)

Parameter	Value
Product Name	PSMA5 Rabbit pAb
Remarks/Alias	PSC5; ZETA; PSMA5
Species	Human
GeneID (Human)	5686
GeneID	5686
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human PSMA5 (NP_002781.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28066
Protein Weight	26kDa
Shipping	Ice bag