

Anti-Psma5 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P28066

Item Number: CM0018017



Introduction

Rabbit polyclonal antibody targeting human PSMA5 (Proteasome subunit alpha type-5), validated for WB, IP, and ELISA. Reacts with human, mouse, and rat. Suitable for proteasome and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	PSMA5 Rabbit pAb
Remarks/Alias	PSC5; ZETA; PSMA5
Species	Human
Gene ID (Human)	5686
Gene ID	5686
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human PSMA5 (NP_002781.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28066
Protein Weight	26kDa
Shipping	Ice bag