


Anti-Psma6 Polyclonal Antibody For Wb, If/Icc, Elisa - P60900

Item Number: CM0021261



Introduction

Rabbit polyclonal antibody targeting human PSMA6 (proteasome subunit alpha type-6), validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for proteasome and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	PSMA6 Rabbit pAb
Remarks/Alias	IOTA; p27K; PROS27; PSMA6
Species	Human
Gene ID (Human)	5687
Gene ID	5687
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human PSMA6 (NP_002782.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60900
Protein Weight	27kDa
Shipping	Ice bag