


Anti-Psma6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P60900

Item Number: CM0014257



Introduction

Rabbit polyclonal antibody targeting human PSMA6, validated for western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Suitable for proteasome and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	PSMA6 Rabbit pAb
Remarks/Alias	IOTA; p27K; PROS27; PSMA6
Species	Human
GeneID (human)	5687
GeneID	5687
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human PSMA6 (NP_002782.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60900
Protein Weight	27kDa
Shipping	Ice bag