


Anti-Psma7 Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, If/Icc, Elisa - O14818

Item Number: CM0006566



Introduction

Rabbit monoclonal antibody against human PSMA7 (Proteasome subunit alpha type-7), KD validated. Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for proteasome and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] PSMA7 Rabbit mAb
Remarks/Alias	C6; HSPC; RC6-1; XAPC7; HEL-S-276; [KD Validated] PSMA7
Species	Human
Gene ID (Human)	5688
Gene ID	5688
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-248 of human PSMA7 (NP_002783.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14818
Protein Weight	28kDa
Shipping	Ice bag