

Anti-Psma/Folh1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q04609

Item Number: CM0007574



Introduction

Rabbit monoclonal antibody against human PSMA/FOLH1. Validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Target is a folate hydrolase and NAALADase involved in prostate cancer and neurotransmission.

[Learn More](#)

Parameter	Value
Product Name	PSMA/FOLH1 Rabbit mAb
Remarks/Alias	PSM; FGCP; FOLH; GCP2; PSMA; mGCP; GCPII; NAALAD1; PSMA/FOLH1
Species	Human
GeneID (Human)	2346
GeneID	2346
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human PSMA/FOLH1 (NP_004467.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04609
Protein Weight	84kDa
Shipping	Ice bag