

Anti-Psma/Folh1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q04609

Item Number: CM0009797



Introduction

PSMA/FOLH1 rabbit monoclonal antibody for Western blot, immunofluorescence, and ELISA. Targets human and rat glutamate carboxypeptidase 2 (Q04609), a key enzyme in folate uptake and prostate tumor progression.

[Learn More](#)

Parameter	Value
Product Name	PSMA/FOLH1 Rabbit mAb
Remarks/Alias	PSM; FGCP; FOLH; GCP2; PSMA; mGCP; GCPII; NAALAD1; PSMA/FOLH1
Species	Human
GeneID (human)	2346
GeneID	2346
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human PSMA/FOLH1 (Q04609).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q04609
Protein Weight	84kDa
Shipping	Ice bag