

Anti-Psmb1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P20618

Item Number: CM0018104



Introduction

PSMB1 Rabbit Polyclonal Antibody (HC5/PSC5) for Western blot, IF/ICC, ELISA. Detects human, mouse, rat proteasome subunit beta type-1 (26kDa). Ideal for proteasome research and protein degradation studies.

[Learn More](#)

Parameter	Value
Product Name	PSMB1 Rabbit pAb
Remarks/Alias	HC5; PSC5; PMSB1; NEDMHAL; PSMB1
Species	Human
GeneID (Human)	5689
GeneID	5689
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-241 of human PSMB1 (NP_002784.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20618
Protein Weight	26kDa
Shipping	Ice bag