


Anti-Psmb10 Rabbit Monoclonal Antibody For Wb, Elisa - P40306

Item Number: CM0003033



Introduction

PSMB10 Rabbit mAb is a specific monoclonal antibody for detecting proteasome subunit beta-2i in human, mouse, and rat samples via WB and ELISA. Suitable for immunoproteasome research, antigen processing, and disease studies.

[Learn More](#)

Parameter	Value
Product Name	PSMB10 Rabbit mAb
Remarks/Alias	LMP10; MECL1; PRAAS5; beta2i; PSMB10
Species	Human
GeneID (Human)	5699
GeneID	5699
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human PSMB10 (P40306).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40306
Protein Weight	29 kDa
Shipping	Ice bag