


Anti-Psmb10 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P40306

Item Number: CM0030094



Introduction

Rabbit polyclonal antibody targeting human PSMB10 (LMP10, MECL1, beta2i), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Suitable for immunoproteasome and antigen processing research.

[Learn More](#)

Parameter	Value
Product Name	PSMB10 Rabbit pAb
Remarks/Alias	LMP10; MECL1; PRAAS5; beta2i; PSMB10
Species	Human
GeneID (Human)	5699
GeneID	5699
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human PSMB10 (NP_002792.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40306
Protein Weight	29kDa
Shipping	Ice bag