

Anti-Psmb2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P49721

Item Number: CM0018511



Introduction

PSMB2 Rabbit pAb is a polyclonal antibody targeting human proteasome subunit beta type-2, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with rat. Ideal for proteasome research and protein degradation studies.

[Learn More](#)

Parameter	Value
Product Name	PSMB2 Rabbit pAb
Remarks/Alias	HC7-I; PSMB2
Species	Human
GeneID (Human)	5690
GeneID	5690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human PSMB2 (NP_002785.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	P49721
Protein Weight	23kDa
Shipping	Ice bag