


Anti-Psmb2 Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P49721

Item Number: CM0030578



Introduction

Rabbit polyclonal antibody targeting PSMB2, a 20S proteasome subunit involved in ubiquitin-dependent protein degradation. Validated for WB, IHC-P, IP, ELISA. Cross-reacts with human, mouse, rat. Essential for proteasome research in cancer, neurodegeneration, and immunology.

[Learn More](#)

Parameter	Value
Product Name	PSMB2 Rabbit pAb
Remarks/Alias	HC7-I; PSMB2
Species	Human
Gene ID	5690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human PSMB2 (NP_002785.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49721
Protein Weight	23kDa
Shipping	Ice bag