


## Anti-Psmb3 Polyclonal Antibody For Wb, Ihc-P, Elisa - P49720

Item Number: CM0025399



### Introduction

PSMB3 Rabbit polyclonal antibody (HC10-II), suitable for Western blot, immunohistochemistry, and ELISA. Detects human, mouse, and rat PSMB3. Essential for proteasome research.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PSMB3 Rabbit pAb                                                                                                                        |
| Remarks/Alias    | HC10-II; PSMB3                                                                                                                          |
| Species          | Human                                                                                                                                   |
| Gene ID (Human)  | 5691                                                                                                                                    |
| Gene ID          | 5691                                                                                                                                    |
| Immunogen        | Recombinant protein   Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human PSMB3 (NP_002786.2). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, IHC-P, ELISA                                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                       |
| SWISS            | P49720                                                                                                                                  |
| Protein Weight   | 23kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |