

Anti-Psmb4 Rabbit Polyclonal Antibody For Wb And Elisa - P28070

Item Number: CM0027319



Introduction

CamelBio supplies a rabbit polyclonal antibody against human PSMB4 (proteasome subunit beta type-4), validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying protein degradation and ubiquitin-proteasome pathways.

[Learn More](#)

Parameter	Value
Product Name	PSMB4 Rabbit pAb
Remarks/Alias	HN3; HsN3; PRAAS3; PROS26; PROS-26; PSMB4
Species	Human
Gene ID	5692
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human PSMB4 (NP_002787.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28070
Protein Weight	29kDa
Shipping	Ice bag