


Anti-Psmb5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P28074

Item Number: CM0022143



Introduction

Rabbit polyclonal antibody to human PSMB5 (proteasome subunit beta type-5). Applications: WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Suitable for proteasome and protein degradation studies.

[Learn More](#)

Parameter	Value
Product Name	PSMB5 Rabbit pAb
Remarks/Alias	MB1, LMPX, PSMB5
Species	Human
Gene ID (Human)	5693
Gene ID	5693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-263 of human PSMB5 (NP_002788.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28074
Protein Weight	28kDa
Shipping	Ice bag