

## Anti-Psmb6 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28072

Item Number: CM0028995



### Introduction

Rabbit polyclonal antibody targeting human PSMB6, validated for Western blot, IHC-P, IF/ICC, and ELISA. Detects PSMB6 in human, mouse, and rat. Suitable for proteasome and protein degradation studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | PSMB6 Rabbit pAb                                                                                                  |
| Remarks/Alias    | LMPY; DELTA; PSMB6                                                                                                |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 5694                                                                                                              |
| GeneID           | 5694                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 35-239 of human PSMB6 (NP_002789.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P28072                                                                                                            |
| Protein Weight   | 25kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |