


## Anti-Psmb7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99436

Item Number: CM0014381



### Introduction

Rabbit polyclonal antibody against human PSMB7 (Proteasome subunit beta type-7), validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Suitable for human, mouse, and rat samples.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PSMB7 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | PSMB7; Z                                                                                                                               |
| Species          | Human                                                                                                                                  |
| Gene ID (Human)  | 5695                                                                                                                                   |
| Gene ID          | 5695                                                                                                                                   |
| Immunogen        | Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 44-200 of human PSMB7 (NP_002790.1). |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, IF/ICC, ELISA                                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                      |
| SWISS            | Q99436                                                                                                                                 |
| Protein Weight   | 30kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |