


Anti-Psmb8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P28062

Item Number: CM0028302



Introduction

Rabbit polyclonal antibody against human PSMB8, a component of the immunoproteasome involved in antigen processing. Validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Suitable for apoptosis and immune response research.

[Learn More](#)

Parameter	Value
Product Name	PSMB8 Rabbit pAb
Remarks/Alias	JMP; ALDD; LMP7; NKJO; D6S216; PRAAS1; PSMB5i; RING10; D6S216E; PSMB8
Species	Human
GeneID (Human)	5696
GeneID	5696
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 163-272 of human PSMB8 (NP_004150.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28062
Protein Weight	30kDa
Shipping	Ice bag