


Anti-Psmb9/Lmp2 Monoclonal Antibody For Wb, Ihc-P, Elisa - P28065

Item Number: CM0011562



Introduction

Rabbit monoclonal antibody against human PSMB9/LMP2 immunoproteasome subunit, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Suitable for antigen processing and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	PSMB9/LMP2 Rabbit mAb
Remarks/Alias	LMP2; PRAAS3; PSMB6i; RING12; beta1i; PSMB9/LMP2
Species	Human
Gene ID (Human)	5698
Gene ID	5698
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28065
Protein Weight	23kDa
Shipping	Ice bag