


## Anti-Psmb9/Lmp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P28065

Item Number: CM0021262



### Introduction

Rabbit polyclonal antibody against human PSMB9/LMP2 (Proteasome subunit beta type-9). Validated for WB, IF/ICC, ELISA. Suitable for immunoproteasome and antigen processing research.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | PSMB9/LMP2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | LMP2; PRAAS3; PSMB6i; RING12; beta1i; PSMB9/LMP2                                                                      |
| Species          | Human                                                                                                                 |
| Gene ID          | 5698                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-219 of human PSMB9/LMP2 (NP_002791.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                                     |
| Cross Reactivity | Human                                                                                                                 |
| SWISS            | P28065                                                                                                                |
| Protein Weight   | 23kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |