

Anti-Psmc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P62191

Item Number: CM0014918



Introduction

PSMC1 Rabbit pAb detects human, mouse, rat PSMC1 in WB, IHC-P, ELISA. Ideal for 26S proteasome and ubiquitin-proteasome pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PSMC1 Rabbit pAb
Remarks/Alias	S4; p56; RPT2; P26S4; NEDGTH; PSMC1
Species	Human
Gene ID (Human)	5700
Gene ID	5700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-440 of human PSMC1 (NP_002793.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62191
Protein Weight	49kDa
Shipping	Ice bag