

Anti-Psmc2 Rabbit Polyclonal Antibody For Wb, Elisa - P35998

Item Number: CM0020582



Introduction

High-quality anti-PSMC2 rabbit polyclonal antibody validated for WB and ELISA; detects human, mouse, and rat 26S proteasome regulatory subunit 7 (S7/MSS1). Ideal for ubiquitin-proteasome pathway and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	PSMC2 Rabbit pAb
Remarks/Alias	S7; MSS1; RPT1; Nbla10058; PSMC2
Species	Human
Gene ID (Human)	5701
Gene ID	5701
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human PSMC2 (NP_002794.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35998
Protein Weight	49kDa
Shipping	Ice bag