


Anti-Psmc3 Rabbit Polyclonal Antibody For Wb - P17980

Item Number: CM0013765



Introduction

PSMC3 Rabbit pAb is a rabbit polyclonal antibody validated for WB, IF/ICC, IP, and ELISA. It targets human PSMC3, a 26S proteasome regulatory subunit. Ideal for proteasome and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	PSMC3 Rabbit pAb
Remarks/Alias	RPT5; TBP1; DCIDP; PSMC3
Species	Human
Gene ID (Human)	5702
Gene ID	5702
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human PSMC3 (NP_002799.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P17980
Protein Weight	49kDa
Shipping	Ice bag