


Anti-Psmc4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P43686

Item Number: CM0007459



Introduction

Validated anti-PSMC4 rabbit monoclonal antibody for WB, IHC-P, ELISA. Detects human, mouse, rat PSMC4, a 26S proteasome regulatory subunit involved in ATP-dependent protein degradation. Ideal for proteasome research.

[Learn More](#)

Parameter	Value
Product Name	PSMC4 Rabbit mAb
Remarks/Alias	S6; RPT3; TBP7; TBP-7; MIP224; PSMC4
Species	Human
GeneID (Human)	5704
GeneID	5704
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human PSMC4 (NP_006494.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43686
Protein Weight	47kDa
Shipping	Ice bag