


Anti-Psmc4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P43686

Item Number: CM0019060



Introduction

Rabbit polyclonal antibody against human PSMC4 (S6/RPT3/TBP7), suitable for WB, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Ideal for researching 26S proteasome regulatory subunit 6B functions in ubiquitin-dependent protein degradation.

[Learn More](#)

Parameter	Value
Product Name	PSMC4 Rabbit pAb
Remarks/Alias	S6, RPT3, TBP7, TBP-7, MIP224, PSMC4
Species	Human
GeneID	5704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human PSMC4 (NP_006494.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43686
Protein Weight	47kDa
Shipping	Ice bag