


Anti-Psmc5 Monoclonal Antibody For Wb, If/Icc, Elisa - P62195

Item Number: CM0003007



Introduction

Rabbit monoclonal antibody targeting PSMC5 (26S proteasome regulatory subunit 8) for WB, IF/ICC, ELISA. Detects human, mouse, rat PSMC5. Ideal for studying the ubiquitin-proteasome system.

[Learn More](#)

Parameter	Value
Product Name	PSMC5 Rabbit mAb
Remarks/Alias	S8; p45; RPT6; SUG1; SUG-1; TBP10; TRIP1; p45/SUG; PSMC5
Species	Human
GeneID (Human)	5705
GeneID	5705
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62195
Protein Weight	46kDa
Shipping	Ice bag