

Anti-Psmc5 Rabbit Polyclonal Antibody For Wb, Elisa - P62195

Item Number: CM0013120



Introduction

Polyclonal antibody targeting human PSMC5 (26S proteasome regulatory subunit 8). Validated for WB and ELISA; cross-reacts with mouse. Ideal for studying proteasome function, protein degradation, and related pathways.

[Learn More](#)

Parameter	Description
Product Name	PSMC5 Rabbit pAb
Remarks/Alias	S8; p45; RPT6; SUG1; SUG-1; TBP10; TRIP1; p45/SUG; PSMC5
Species	Human
GeneID (Human)	5705
GeneID	5705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-406 of human PSMC5 (NP_002796.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P62195
Protein Weight	46kDa
Shipping	Ice bag