


Anti-Psmc5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62195

Item Number: CM0014551



Introduction

High-quality rabbit polyclonal antibody targeting human PSMC5, a 26S proteasome regulatory subunit. Validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for proteasome and ubiquitin-proteasome pathway research.

[Learn More](#)

Parameter	Value
Product Name	PSMC5 Rabbit pAb
Remarks/Alias	S8; p45; RPT6; SUG1; SUG-1; TBP10; TRIP1; p45/SUG; PSMC5
Species	Human
GeneID (Human)	5705
GeneID	5705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-406 of human PSMC5 (NP_002796.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62195
Protein Weight	46kDa
Shipping	Ice bag