


Anti-Psmc6 Rabbit Polyclonal Antibody For Wb And Elisa - P62333

Item Number: CM0027183



Introduction

High-quality rabbit polyclonal antibody against human PSMC6 (p42/RPT5/SUG2), validated for WB and ELISA. Recognizes human and mouse 26S proteasome regulatory subunit 10B. Ideal for ubiquitin-proteasome pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PSMC6 Rabbit pAb
Remarks/Alias	p42; RPT5; SUG2; PSMC6
Species	Human
GeneID (Human)	5706
GeneID	5706
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-150 of human PSMC6 (NP_002797.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P62333
Protein Weight	44kDa
Shipping	Ice bag