

Anti-Psmd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q99460

Item Number: CM0018910



Introduction

Rabbit polyclonal antibody against human PSMD1 (26S proteasome non-ATPase regulatory subunit 1). Validated for Western blot and ELISA applications in human, mouse, and rat samples. Suitable for proteasome and ubiquitin pathway research.

[Learn More](#)

Parameter	Details
Product Name	PSMD1 Rabbit pAb
Remarks/Alias	S1; P112; Rpn2; PSMD1
Species	Human
GeneID (Human)	5707
GeneID	5707
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 794-953 of human PSMD1 (NP_002798.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99460
Protein Weight	106kDa
Shipping	Ice bag